SCHARIOTH
MACULA LENS
A45SML

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1. WHO causes of blindness. 2008. 2. Agarwal A, Lipshitz I, Soosan J, Mandeep L, Rahul T, Kumar AD, Agarwal A.: **Mirror telescopic intraocular lens for age-related macular degeneration.** Design and preliminary clinical results of the Lipshitz macular implant. *J Catar Refr Surg*. 2008; 34: 87–94. 3. A.Calabrèse, JB Bernard, L. Hoffart, G. Faure, F. Barouch, J.Conrath, E. Caster: **Wet versus Dry Age-Related Macular Degeneration in Patients with Central Field Loss:** Different Effects on Maximum Reading Speed. *Investigative Ophthalmol & Visual Science*, April 2011, Vol.52, No. 5. 4. Lipshitz, I, Loewenstein, A, Reingewirtz, M, and Lazar, M. **An intraocular telescopic lens for macular an intraocular telescopic lens for macular degeneration** *Ophthalmic Surg Lasers*. 1977; 28: 513–517. 5. Age-Related Eye Disease Study Research Group. **Potential public health impact of age-related eye disease study results; AREDS report no. 11** *Arch Ophthalmol*. 2003; 121: 1621–1624. 6. Rosner, M., Ben-Simon, G., and Sachs, D. **Feasibility and safety of laser treatments in eyes with an intraocular implantable miniature telescope.** *J Cataract Refract Surg*. 2003; 29: 1005–1010. 7. Agarwal, A. **Consultation section: cataract surgical problems.** *J Cataract Refract Surg*. 2007; 33: 183–185. 8. Riusala, A., Sarna, S., and Immonen, I. **Visual Function Index (VF-14) in exudative age-related macular degeneration of long duration.** *Am J Ophthalmol*. 2003; 135: 206–212. 9. Kaskaloglu, M., Üretmen, Ö, and Yağcı, A. **Medium-term results of implantable miniaturized telescopes in eyes with age-related macular degeneration.** *J Cataract Refract Surg*. 2001; 27: 1751–1755. 10. Alió, J.L., Mulet, E.M., Ruiz-Moreno, J.M. et al. **Intraocular telescopic lens evaluation in patients with age-related macular degeneration.** *J Cataract Refract Surg*. 2004; 30: 1177–1189. 11. Koziol, J, Peyman, G.A, Cionni, R et al. **Evaluation and implantation of a teledioptric lens system for cataract and age-related macular degeneration.** *Ophthalmic Surg*. 1994; 25: 675–684. 12. Peli, E. **The optical functional advantages of an intraocular low-vision telescope.** *Optom Vis Sci*. 2002; 79:225–233. 13. Goldstein, H. **The reported demography and the causes of blindness throughout the world.** *Adv Ophthalmol*. 1980; 40: 1–99. 14. Hengerer FH, Artal P, Kohnen T, Conrad-Hengerer I. **Initial Clinical Results of a New Telescopic IOL Implanted in Patients With Dry Age-Related Macular Degeneration.** *J Refract Surg*. 2015;31(3):158–162. 15. K. Michalska-Malecka, A. Kabiesz, M. Nowak, D. Spiewak: **Age related macular degeneration – challenge for future:** Pathogenesis and new perspectives for treatment. 16. Friedman DS, O’Colmain BJ, Muñoz B, Tomany SC, McCarty C, de Jong PT, Nemesure B, Mitchell P, Kempen J. Eye Diseases Prevalence Research Group. **Prevalence of age-related macular degeneration in the United States.** *Arch Ophthalmol*. 2004;122:564–572. 17. Vingerling JR, Dielemans I, Hofman A, Grobbee DE, Hijmering M, Kramer CF, de Jong PT. **The prevalence of age-related maculopathy in the Rotterdam Study.** *Ophthalmology*. 1995;102:205–210. 18. Klein R, Knudtson MD, Lee KE, Gangnon RE, Klein BE. **Age-period-cohort effect on the incidence of age-related macular degeneration: the Beaver Dam Eye Study.** *Ophthalmology*. 2008;115:1460–1467. 19. Klein R, Klein BE, Lee KE, Cruickshanks KJ, Gangnon RE. **Changes in visual acuity in a population over a 15-year period: the Beaver Dam Eye Study.** *Am J Ophthalmol*. 2006;142:539–549. 20. Parmeggiani F, Romano MR, Costagliola C, Semeraro F, Incorvaia C, D’Angelo S, Perri P, De Palma P, De Nadai K, Sebastiani A. **Mechanism of inflammation in age-related macular degeneration.** *Mediators Inflamm*. 2012;2012:546786. 21. Tien Y. Wong, Xiang Li, Xinyi Su, Gemmy C. Cheung, Ching-Yu Cheng: **The number and distribution of People with Age-Related Macular degeneration Worldwide in 2020 and 2040:** A Systematic Review and Hierarchical Bayesian Meta-Analysis. Program NO: 220; Poster Board Nu:DOO65; ARVO 2013. 22. Wan Ling Wong, Xinyi Su, Xiang Li, Chui Ming G Cheung, Ronald Klein, Ching-Yu Cheng, Tien Y. Wong: **Global prevalence of age related macular degeneration and disease burden projection for 2020 an 2040: a systematic review and meta-analysis.** *Lancet Glob Health* 2014;2:e106–16; publishe on-line January 3, 2014. 23. Yalin Mu, Manli Zhao, Guangming Su: **Stem cell-based therapies for age-related macular degeneration: current status and prospects.** *Int J Clin Exp Med* 2014;7(11):3843–3852. 24. Scharioth GB, “Add-on IOLs for Vitrectomized Silicone oil-filled Eyes,” *Retinal Physician* May 2014; 11:31–33. Available at: <http://www.retinalphysician.com/articleviewer.aspx?articleID=111129>. 25. Sauder G, Cordes A. *Scope of applications and experience record with a new generation of add-on IOLs: the A4W lens.* *Ophthalmo-Chirurgie* 2012; 24(suppl. 2):1–6. Available at: http://www.1stq.de/media/raw/A4W_Englisch_AddOn_OphthalmoChirurgie_engl.pdf. 26. Scharioth GB (2015): **New add-on intraocular lens for patients with age-related macular degeneration.** *J Cataract Refract Surg* (in Print).