

# LESS STEPS – MORE VISION



## **TRIFOCAL**

Complete Spectacle Independence

## **7 RINGS TO PERFECTION**

Superior Contrast Sensitivity

**LIBERTY**®

**MEDICENTUR**

Material. Design. Optics.



[liberty-lens.com](http://liberty-lens.com)



## ENHANCED QUALITY OF VISION

IOL manufacturers eagerly compete for superior intermediate performance or compress the addition towards the intermediate to create some „comfort“. While doing this they disregard the main key of success for multifocal IOLs – SPECTACLE INDEPENDENCE.

**Medicontur chooses not to compromise IOL performance.**

Liberty equates to spectacle independence without the compromise of reduced image quality or visual acuity at near and far distances.

**LIBERTY  
FOR YOUR PATIENTS**



## Less Diffractive Steps – More Vision

Light is scattered within an IOL by every diffractive step. Considerable additional light scattering is caused by the manufactured imperfections of these steps on multiple points.

Therefore **not only the quality but also the quantity of the manufactured steps has a great impact on light scattering and loss of light energy.**

The below symbolic figures demonstrate the difference in light scattering and energy loss on two separate diffractive IOL profiles.

### Quantity of Scattered Light




Figure 1a:  
LIBERTY optic with  
**7 diffractive steps**  
Moderate light  
scattering




Figure 1b:  
Diffractive optic with  
**2× more diffractive steps**  
creates minimum  
**3× more light scattering**

What does **LESS STEPS** mean for your patients?

Less light scattering improves vision by


- **Higher Contrast Sensitivity**
- **Less Halos & Glare**



# Out-Of-Focus Light Intensity Distribution by Visualization of Light Propagation along the Optical Axis

Light Energy Intensity at 4.5 mm Aperture\*

→ Halos in FAR focal plane




## Which Lens Would You Prefer For Driving At Night?

\* Simulated by Zemax extension for calculation of light propagation after custom defined diffractive surfaces.




# A TRIFOCAL DESIGN: Elevated Phase Shift\* (EPS) technology

Medicontur's proprietary, patented approach to trifocality uses **EPS** in the central part of the optic to generate constructive interference between the 0th (far) and 1st (near) diffractive order. This **revolutionary optical design** generates the 3rd focal point for intermediate vision in a very unique way:



**Figure 1:** Radial profile with central diffractive phase shift elevated from a destructive level to a constructive interference level.

EPS\* technology uses only **7 rings** within a **3 mm** diameter for complete spectacle independence, improved image quality and excellent performance at near and far.



**Figure 2:** MTF through focus curves with central diffractive phase shifts at a destructive and a constructive interference level resulting in a 3rd focal point for intermediate vision.




# 7 steps to Trifocal perfection

## LIBERTY IOLs achieve trifocal performance with

- **7 steps** in a precise diffractive array within a
- **3 mm** diameter leaving a
- **75%** refractive lens surface.

**LIBERTY IOLs** are strongly **pupil dependent** using the **NEAR TRIAD** which implies miosis under accommodation. We believe that too much light distribution into the near focus above 3 mm aperture does not match normal ocular physiology.

The 4 charts below show the **Useful Light Distribution in % of LIBERTY** and 3 competitors depending on the aperture [mm] \*. **LIBERTY provides the highest light distribution in the near focus under accommodation and the highest for far vision under scotopic conditions.**



\* Based on Strehl ratio calculated from Zemax simulated MTF values

# LIBERTY<sup>®</sup>

"The most balanced trifocal IOL on the market."

M. Assouline, MD, PhD  
France

APAO 2018,  
Hong Kong





**10**  
CLINICAL  
STUDIES

IN  
**12**  
COUNTRIES

MORE THAN  
**900**  
EYES

1. TRIFOCAL PERFORMANCE \*
2. COMPLETE SPECTACLE INDEPENDENCE \*
3. EXCELLENT FAR – NEAR & GOOD INTERMEDIATE VISION \*
4. REFRACTIVE PREDICTABILITY & LONG TERM STABILITY \*
5. OUTSTANDING CONTRAST SENSITIVITY \*
6. EXCELLENT READING SPEED \*
7. MINIMAL REPORTS OF DYSPHOTOPSIA \*



\* Confirmed by studies; see references.

■ Very Satisfied ■ Satisfied ■ Rather satisfied






THE MOST BALANCED TRIFOCAL IOL




# TRIFOCAL performance and Spectacle Independence

Clinical evidence shows excellent far & near, good intermediate vision.

## Outstanding Visual Performance for Liberty 677MY (Study C)




## Equivalent Clinical Visual Acuity compared to FineVision (PhysIOL) (Study B)







# High Image QUALITY


## Superior Contrast Sensitivity compared to FineVision (PhysIOL) in different light conditions. (Study B, C)



## Equivalent Contrast Sensitivity to Bi-Flex Monofocal IOL (Study E)



## Equivalent Reading Speed to Bi-Flex Monofocal IOL (Study E)




„After over 100 implantations trifocal performance is evident. Less halos & glare, better contrast sensitivity. All my patients implanted with Liberty are spectacle independent.“ **J. Györy, MD, Hungary** › ESCRS 2017, Lisabon

# Refractive Stability and Predictability

Excellent Long Term Refractive Stability and Visual Outcomes. (Study C)

## UCVA vs BCVA



Cumulative Snell Visual Acuity (20/40 or better)

## ACCURACY & PREDICTABILITY

Mean spherical equivalent residual error of  $-0.15D (\pm 0.33D)$  (Based on studies A, C, D)



Accuracy of refractive outcome of Liberty overshadow distinctively literature average =  $\pm 0.5 D$ : 75.1%.

Cooke DL, Cooke TL. J Cataract Refract Surg. 2016;42(8):1157-1164



# Liberty 677MY / 677PMY Technical Specification



<b>Type</b>	Single-piece hydrophilic intraocular lens for implantation into the capsular bag
<b>Material</b>	25% water content with UV absorber + blue light filter; Refractive Index 1.46; ABBE number 58
<b>Optic design</b>	Biconvex, aspheric, diffractive-refractive, apodized
<b>Powers available*</b>	+8.00 D → +35.00 D (increment: 0.50 D)
<b>Diffractive zone</b>	EPS technology**, anterior surface, Ø 3.0 mm
<b>Addition</b>	+3.50 D (near); +1.75 D (intermediate)
<b>A-constant***</b>	118.9 (SRK/T)
<b>Dimensions</b>	Overall length 13.0 mm; optic Ø 6.0 mm
<b>PCO protection</b>	360° Special Square Edge (patented)
<b>Haptic angulation</b>	0° – asymmetrical design with posterior vaulting
<b>Sterilization</b>	Steam (shelf life 5 years after sterilization)
<b>Storage conditions</b>	at +15 °C – +35 °C (15% – 50% humidity)

\* Other powers upon request

\*\* Patent pending

\*\*\* It is recommended that surgeons personalize the constants they use.

## Single use injection system for Liberty 677PMY preloaded IOLs

### 1<sup>st</sup> CLICK LOADING


1. Put a drop of visco in the lens holder
2. Insert the injector into the wet container



### 2<sup>nd</sup> CLICK CLOSING


3. Inject viscoelastic material into the cartridge nozzle

4. Fold up the cartridge nozzle




### 3<sup>rd</sup> CLICK READY

5. Push the red stopper forward to the end position
6. Fill visco



### INJECTION

7. Push the plunger midway then pull back to the end and continue pushing



# Material. Design. Optics.

## PREMIUM MATERIAL & ADVANCED MANUFACTURING TECHNOLOGY

## A LEADING HAPTIC DESIGN FOR ROTATIONAL STABILITY

### SAFETY RECORDS

- 30 years experience in IOL design and manufacturing

### IMAGE QUALITY

- Highest ABBE no. (58) for the lowest chromatic aberration

### PCO protection

- Low ionicity surface for less cell adhesion

### RETINA PROTECTION


- UV blocker
- Violet light filtering

### PCO PROTECTION

- 360° posterior optic edge ( $\leq 10$   $\mu\text{m}$  edge radius)

### ROTATIONAL STABILITY

- Ergonomic-adaptive fit with optimal balance in haptic force – DOUBLE LOOP HAPTIC
- 180° contact angle between haptic and the capsular bag equator at 9 mm  $\varnothing$



## TRIFOCAL PERFORMANCE – SPECTACLE INDEPENDENCE

### HIGH QUALITY OF VISION

- 7 steps in a precise diffractive array with a
  - 3 mm diameter leaving a
  - 75% refractive lens surface
- Minimal dysphotopic phenomena
  - Superior contrast sensitivity
  - Uncompromised near vision
- PUPIL DEPENDENT using physiological NEAR TRIAD

**LIBERTY**<sup>®</sup>

**MEDICENTUR**

Material. Design. Optics.

# Liberty 677MTY Technical Specification



<b>Type</b>	Single-piece hydrophilic trifocal toric intraocular lens for implantation into the capsular bag
<b>Material</b>	25% water content with UV absorber + blue light filter; Refractive Index 1.46; ABBE number 58
<b>Optic design</b>	Biconvex, aspheric, diffractive-refractive, apodized, toric
<b>Powers available*</b>	+8.00 D → +35.00 D (increment: 0.50 D)
<b>Cylinders available**</b>	+1.00 D → +4.50 D (increment: 0.50 D) +5.25 D; +6.00 D
<b>Diffractive zone</b>	EPS technology***, anterior surface, Ø 3.0 mm
<b>Addition</b>	+3.50 D (near); +1.75 D (intermediate)
<b>A-constant****</b>	118.9 (SRK/T)
<b>Dimensions</b>	Overall length 13.0 mm; optic Ø 6.0 mm
<b>PCO protection</b>	360° Special Square Edge (patented)
<b>Haptic angulation</b>	0° – asymmetrical design with posterior vaulting
<b>Sterilization</b>	Steam (shelf life 5 years after sterilization)
<b>Storage conditions</b>	at +15 °C – +35 °C (15% – 50% humidity)

\* Other powers upon request

\*\* Please note that at SEQ.: +8.0 D to +9.5 D, only cyl. +1.0 D to +4.5 D applies.

\*\*\* Patent pending

\*\*\*\* It is recommended that surgeons personalize the constants they use.

**LIBERTY**<sup>®</sup>

**MEDICENTUR**


## References

- A.** AF, Dunai et al. (Hungary) *Comparison of two multifocal IOL Types (50 Bi-Flex M & 50 Acrysof IQ) – long term visual outcome.* ESCRS 2016, 2017. **B.** E. Van Acker, MD, (Belgium) *Comparison of clinical outcomes and patient satisfaction after implantation of two different types of diffractive apodized IOLs: 40 Bi-Flex M (based on PAD technology) 40 Finevision (Micro F) trifocal diffractive IOL. Prospective, randomized, observational study.* ESCRS 2017. **C.** J. Györy (Hungary) *Long term functional and morphological outcomes and patient satisfaction after cataract surgery with 100 BiFlex M implantation with/without posterior central circular capsulorhexis (PCCC).* ESCRS 2016, 2017. **D.** J. Fernandez MD, PhD, (Spain) *Visual performance of patients implanted with progressive PAD Bi-Flex M analyzed by the Qvision iPAD Multifocal LensAnalyzer.* ESCRS 2017. **E.** E. Law et al. (UK) *Randomised clinical trial of the Bi-Flex M multifocal intraocular lens.* ESCRS 2017. **F.** A Dextl. (Austria) *Visual Outcome, Patient Satisfaction and Spectacle Independency after Implantation of 50 eyes ith Progressive Bi-Flex M. Final Result of a Multicentric Trial with 50 Consecutive Patients.* ESCRS 2015. **G.** C. Naval. (Philippines) *Evolution of multifocal practice.* ESCRS 2014. **H.** N. Meijide (Argentina) *Visual Outcomes After Bilateral Implantation of a New Progressive Apodized Diffractive Trifocal IOL.* Clinical outcome ASCRS 2018 ID 42517. **I.** J. García-Bella et al. (Spain) *Visual outcomes after progressive apodized diffractive intraocular lens implantation.* Eur J Ophthalmol. 2017 Sep. **J.** M. Assouline MD, PhD, (France) *Comparative Outcome of Four Multifocal Intraocular Lens.* ESCRS 2015.

# FOCUS on the patients vision since 1989

Medicontur stands  
for consistent high quality,  
proven by more than

**6 million  
intraocular lenses  
implanted.**



**MEDICONTUR**

**Medicontur**

Medical Engineering Ltd  
export@medicontur.com  
www.medicontur.com

**Head Office**

Herceghalmi Road 1.  
2072 Zsámbék  
Hungary

**International Office**

Chemin des Aulx 18  
1228 Plan-les-Ouates  
Geneva/Switzerland

CE  
1639